

Legend 2 Exigo



The Handicare Legend 2 Exigo is a very well-equipped crossframe wheelchair.

Its intuitive and simple adjustments, superb sitting comfort and high quality ensure that Legend 2 Exigo satisfies requirements for both ease of maintenance and re-use.

Most of the wide range of accessories available are compatible with the rest of Handicare's wheelchair range.



handicare

LEGEND 2 Exigo



Stability



A stable frame is important to achieve good rolling characteristics, both when turning and when rolling forwards. The unique frame construction ensures maximum stability.

Transfer



Easy transfer is a key factor for a successful wheelchair. The Exigo is equipped with flip-up, detachable armrests as well as leg supports which can be swung out and removed.

User stability and comfortable seating



Achieving the correct knee angle is an important factor in keeping the user stable and comfortably seated. Therefore, the leg support hanger is adjustable to 70, 80 and 90 degrees.

Individual fittings

In order to meet the needs of each user the Legend 2 Exigo can be adjusted in numerous ways in its standard configuration, and a wide range of accessories is also available.

Value for money

Many adjustments possible in standard configuration, excellent design and improved quality.

Wide range of accessories

- Angle adjustable legrest
- Angle adjustable back
- Velcro back
- Calf cushion for angle adjustable legrest
- Height adjustable push handles
- Universal table
- Anti tip stabiliser
- Basic neck rest (requires pushing bar)
- Cover for neck rest
- Pushing bar
- Amputee support
- Armrests for hemiplegia
- Assistant brake
- Plaster-cast support
- Crutch holder
- Spoke guards

Seat width	36, 39, 42, 45, 48 and 51 cm
Seat depth	39 cm (adjustable 39 - 44 cm)
Seat height front	36 - 48 cm
Back height	42 - 45 cm
Armrest height	22 - 27 cm
Total height	89 cm
Total width	Seat width + 18 cm
Transport width	24,5 cm
Total length	100 cm
Weight	15,9 kg
Max user weight	140 kg