

 **kuschall**[®]
designed for life



 **kuschall**[®] Ultra-Light

*“There always comes a time when you must
choose between contemplation and action.”*

Albert Camus

küschall® Ultra-Light



The best driving feel

in a crossfoldable wheelchair



Legendary for its *driving performance* and *stability*



Experience the *lightness* combined with a *high-tech design*

*Markus Freda,
T1/T2 incomplete since 1995, Military injury*

"My work, my family and my passion for tennis certainly keep me on the move. I would never want to give up on an active life. My küschall Ultra-Light gives me the possibility to pursue my interests. Just a few quick moves and it's in the car and I'm ready to go."



Options



One piece flip up footplate for more stability

Front castor supporter with two horizontal positions

Foldable push handles

Pushrims: Aluminium (1), Titanium (2), Stainless steel (3), Supergripp (4), Tetra 2000 (5) and Maxgrepp (6)

High-Performance wheel for lightness.

Küschall wheel bag: For two 24", 25" or 26" wheels

These are only a few of the many options available for this model. Contact your nearest küschall specialist for more information.


Technical Data

	Seat width: 32 - 50 cm (in increments of 2 cm)		Seat depth: 32 - 50 cm (in increments of 2 cm)		Seat height: front: 43 - 55 cm rear: 39 - 50 cm (in increments of 1 cm)
	Backrest height: 30 - 51 cm (in increments of 1,5 cm)		Knee-to-heel length: 20 - 23 cm (high mounted) 30 - 52 cm (in increments of 1 cm)		Seat angle: variable, standard 6 cm
	Backrest angle: -8°, -4°, 0°, 4°, 8°, 12°		Total width: Seat width + 17 cm		Total length (standard): 70° frame approx. 89 cm 80° frame approx. 86 cm
	Total width, folded: approx. 28 cm		Total weight/unloaded: approx. 9,8 kg (seat width 36 cm)		Max. user weight: 130 kg (Double-crossbrace for seat widths of 42 cm and more!)
	Transport weight (without rear wheels): approx. 6,1 kg (seat width 36 cm)				

NOTE: Except where otherwise specified, all measurements and weights indicated are based on the standard of this model in seat width 36cm.

Colours

Frame



CPS (Colour Power Set)



CPS includes Starec hubs with quick release protection, 22" and 24" Starec rear wheel rims, adapter plates, forks, castor housing attachment when option is selected.

Küschall AG | Benkenstrasse 260 | CH-4108 Witterswil | Switzerland | Phone +41 (0) 61 487 70 70 | Fax +41 (0) 61 487 70 79 | www.kuschall.com

United Kingdom:

Invacare Ltd. | South Road | Bridgend Industrial Estate | Bridgend CF31 3PY | Phone +44 (0)1 656 647327 | Fax +44 (0)1 656 649016 | www.invacare.co.uk

Ireland:

Invacare Ltd. | Unit 5, Seatown Business Campus | Seatown Road | Swords, Co. Dublin | Phone +353 (0)1 810 7084 | Fax +353 (0)1 810 7085 | www.invacare.ie

küschall® is a registered trademark. Every effort has been made to ensure that the contents of this publication are fully up-to-date at the time of going to print. As part of its ongoing improvement of products, Küschall AG reserves the right to modify existing models at any time. Any use of this publication, or of parts thereof, as well as any reproduction of images, must have the written consent of Küschall AG.