

küschall[®]
designed for life















Simply extreme and superlight

The K_{SL}

The K^{SL}




Technical Data

	Seat width: 34-42 cm in increments of 2 cm		Backrest angle: -8°, -4°, 0°, 4°, 8°
	Seat depth: 37,5 cm / 40 cm 42,5 cm		Total width (with 3° camber): SB + 23 cm
	Seat height: front: 50 cm (24") front: 51 cm (25") rear: 38 - 48 cm in increments of 1 cm		Total length: 75°: ca. 83 cm 90°: ca. 76 cm
	Backrest height: 27 - 39 cm in increments of 1.5 cm		Total weight unloaded: Approx. 7 kg
	Knee-to-heel length: 40-48 cm in increments of 1 cm		Max. user weight: 100 kg
	Seat angle: variable, standard 7 cm		Transport weight w/out rear wheels, Approx. 4,5 kg


NOTE: Except where otherwise specified, all measurements and weights indicated are based on the standard of this model in seat width 36 cm.

Options



These are only a few of the many options available for this model. Contact your nearest kuschall specialist for more information.

Colours



CPS set includes: High Performance wheel hubs with quick release protection, forks, castor housing caps.

A Masterpiece; Lightness and design in perfect harmony.

The lightest chair on the market: 4.5 kg (without rear wheels)



Freedom through supreme lightness.

**Oliver Anthofer,
Paraplegia, Th 12/L1**

Freedom is key. The KSL by küschall is the culmination of design, performance and extreme lightness giving me the perfect tool for complete independence. With it I can achieve anything I set my mind to.

*Life is a succession
of moments. To live
each one is to
succeed.* Corita Kent