



Breezy 100



BREEZY

**INNOVATE.
BUILD FOR THE FUTURE**



Models

An entry-level steel chair, the Breezy 100 range is an economical choice for a rugged and durable wheelchair. Armrests can be flipped back and removed for ease of transfer. Removable footrests, which can swing in and out for those areas where space is limited

All of the critical functional parts of the chair have been designed to make it much easier for users to operate. The mechanism for the removal of the footrest is large, accessible and very simple for users or carers to operate.

Available in seat widths of 15, 16, 17, 18, 19 or 20" and green, chrome or metallic blue frame colours. The maximum user weight is 19.7 stone.



115



110



Breezy 100 with Folding backrest

The half folding backrest folds down for easy storage and make it easier to fit your wheelchair into the boot of your car.

Breezy 115:
600 mm (24")
solid rear wheel

Breezy 110:
300 mm (12")
solid rear wheel



121



141

Breezy 100 with Reclining backrest

For those looking for extra comfort, the Breezy 121/141 offers the same features as the Breezy 110/115, but with a reclining back for you to relax in when tired..

Breezy 121:
600 mm (24")
solid rear wheel

Breezy 141:
300 mm (12")
solid rear wheel

BREEZY

Breezy 100

Options and accessories



Ref. **12 E**
Long armrest



Ref. **15 E**
Height adjustable armrest



Ref. **74 E**
Belt

Ref. **98 E**

Double amputee



Ref. **21 E**

ELR



Ref. **71 E**

Tray Table



Ref. **72 E**

Crutch support



Ref. **32 E**

Peumatic wheels 24"



Ref. **35 E**

Quick release



Ref. **95 E**

Drum brakes



Ref. **77 E**

Anti tips



Ref. **36 E**

Hemiplegic



Ref. **70 E**

Wheel lock extensions



Ref. **55 E**
Standard headrest



BREEZY

Dimensions:

Seat		Armrests		Backrest		General measurements				Standard
										Total WEIGHT
Model 115										
37,40, 43,46 (*)	51	41	22	45	37,40, 43,46 (*)	55,58, 61,64	94	105	30	17,9
Model 110										
37,40, 43,46 (*)	51	41	22	45	37,40, 43,46 (*)	51,54, 57,60	94	95	28	16,5
Model 121										
37,40, 43,46 (*)	51	41	22	82	37,40, 43,46 (*)	55,58, 61,64	128	115	31	20,6
Model 141										
37,40, 43,46 (*)	51	41	22	82	37,40, 43,46 (*)	51,54, 57,60	128	106	28	18,6

Colours: Chrome-plated, metallic blue, metallic green

Maximum user weight - 125kgs (19.7 stone)

(*) Widths of 49 and 52 cm are available as special sizes

Tolerance: ± 2 cm. (measurements in cm.)

Standard characteristics: Breezy 100

- Removable footrests, which can be swung inwards and outwards
- Armrests, which can be flipped back and removed
- Upholstered in leather feel upholstery
- Threaded axles as standard
- Solid 200 mm (8") front wheel
- Solid 600 mm (24") sport rear wheel/solid 300 mm (12") wheel

STPS
Sunrise Training and Education Programmes

KNOWLEDGE FOR THE FUTURE

**SUNRISE
MEDICAL**

Sunrise Medical Limited
High Street
Wollaston, Stourbridge
West Midlands DV8 4PS, England
Phone +44 (0) 1384 446622
Fax +44 (0) 1384 446644
www.sunrisemedical.com

**INNOVATE.
BUILD FOR THE FUTURE**