



# Emerald

Standard for comfort and drive



# STERLING

INNOVATE.  
BUILD FOR THE FUTURE

# STERLING



# Emerald



Be inspired by the Emerald a scooter for life ...an outdoor life

If you are looking for style & comfort in a compact package, consider the new Sterling Emerald. The Emerald has a sporty look with a lively 6 mph maximum speed. The full suspension chassis offers an exciting ride with great stability. Larger wheels & narrow profile tyres make for easier steering but at the same time easy to drive over rough terrain.



Suspension for additional comfort



Stylish alloy wheels



Stylish lights & indicator package



Fully adjustable captain seat

Overall Length	48" / 122 cm
Overall Width	22" / 55 cm
Maximum user weight	300 lbs / 137 kg / 21 st
Maximum Gradient	14 degrees - 115 kg user weight 11 degrees - 137 kg user weight
Range	22 m / 36 km
Speed	6 mph / 9.5 kph
Turning Radius	55" / 140 cm
Weight (excluding batteries)	158 lbs / 72 kg
Battery Capacity	36 AH
Adjustable seat height	17" - 20" / 43 cm - 51 cm
Front wheel diameter	11" / 28 cm
Rear wheel diameter	11" / 28 cm

Range given represents test conditions. Please note that temperature, user weight, terrain and battery capacity may affect overall performance.

**STLPS**  
Sunrise Training and Education Programmes

KNOWLEDGE FOR THE FUTURE



Sunrise Medical Limited  
High Street Wollaston  
West Midlands DY8 4PS  
England

Phone +44 (0) 1384 446622  
Fax +44 (0) 1384 446644  
[www.sunrisemedical.com](http://www.sunrisemedical.com)



**INNOVATE.  
BUILD FOR THE FUTURE**